

ECGS Quantification by Bradford Method

5/20/16

1.

Dilute 5X Protein Assay Buffer (Bio-Rad, stored at 4C) to 1X

Make enough for X number of samples

e.g. standard 5 samples (BSA high 1:10, med 1:20, low 1:100, ECGS #1, ECGS #2) means

6 mL PAB/ 24 mL ddH₂O (make extra for pipetting errors)

5 mL per sample

2.

Prepare eppendorf tubes

Add 100 uL each sample diluted in H₂O, **add ECGS in hood** (otherwise not sterile)

e.g. for 1:20 dilution add 95 uL H₂O and 5 uL sample

3.

Add each 100 uL prepped sample to 5mL PAB tubes

Vortex immediately

Sit 5 min-1hr

4. Read on plate reader, select absorbance and wells

200 uL sample per well

1X PAB, Sample 1, Sample 2, Sample 3, ... Sample n